49 E / 49 EP / 49 IS

Ref. 6030062-B

Your Endo Specialist™

MICRO-MEGA®

5-12, rue du Tunnel
F-25006 BESANCON Cedex

DISTRIBUTION :

France :
MICRO-MEGA®

Tel : +33 (0)3 81 54 42 15
Fax : +33 (0)3 81 54 42 44
atelier.reparation@micro-mega.com

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I. WARNINGS

1.1 General warnings

- The device is meant for professional dental use only.
- Perform regular maintenance.
- Use the handpiece only for operations of dental drilling. Any other use or alteration of the device is banned and may be dangerous.
- Verify that only instruments in good condition are used.
- Verify that the instrument is properly and firmly held before starting the handpiece.
- Use a motor in good condition.
- Set the speed of the motor in accordance with the utilized handpiece.
- Eliminate diamond drills where grains are partially missing.
- Run the handpiece with instrument in place only.
- Avoid working with a handpiece which heats up abnormally.
- It is recommended to wear appropriate protection gloves and eye-protection while handling the handpiece and the rotary instruments.
- Use a dental dam during intervention on a patient.
- During long treatments, repeat lubrication every ten minutes.
- Contact the MICRO-MEGA® after sales service as soon as a problem with the handpiece occurs.

FR

EN

DE

IT

ES

- Stop using the handpiece in case of damage, overheating, irregular operating sound, excessive vibration, excessive run-out of the instrument or inadequate instrument fixation. In order to prevent further damage, contact the MICRO-MEGA® after sales service or your approved MICRO-MEGA® dealer.




1.2 Warnings for reprocessing (cleaning - disinfection - lubrication - sterilization)

- Sterilize and lubricate the handpiece before first use and after every intervention on a patient.
- Remove the rotary instruments (drills, burs, etc.) before reprocessing.
- Passing the handpiece in a basin or an ultrasound system is banned.
- Soaking the handpiece in a detergent product provided to this effect must not exceed 15 minutes.
- Only use a detergent-disinfectant product specifically intended to the pre-disinfection of dental handpieces and contra-angles which is compatible with aluminium and does not contain soda based products or aldehyde. Do not use any alkaline detergents or solutions (pH > 7).
- Comply with the concentration recommended by the pre-disinfectant manufacturer.
- Respect the immersion time recommended by the pre-disinfectant manufacturer.
- Do not use vinegar to clean the spray nozzles.
- Open the sterilization bag just before use.
- Excessive use of maintenance products (lubricants, cleaning and disinfectant products) may induce the penetration of these products into the inside of the electrical brush motor and affect its functioning. It is imperative to follow the maintenance instructions accompanying each product.
- In case the handpiece is not used for an extended period of time, clean and lubricate it in accordance with the instructions hereafter and store it in a dry place.

II. DISMANTLING AND REASSEMBLY OF THE HANDPIECES 49 E - 49 EP - 49 IS

2.1 Removing the interior from the outer shell

- Insert and lock the supplied pin into the handpiece.
- Slide the wrench without rotating it into the rear as far as it may go. Make sure that the two arrows on the wrench and on the rear of the handpiece are aligned.
- Hold the handpiece vertically and press the pin down against a hard surface in order to unlock the axle.
- Pull the axle from the handpiece shell.

FR

EN

DE

IT

ES

2.2 Removing the twin-chuck set



- Tightly hold forward the spindle by sliding the milled end of the wrench over it and unscrew the rear axle using the flat spanner.
- Remove the twin-chucks and their central ring together with the pin.

2.3 References and reassembling of the twin-chuck set

Ref. 59800078


Twin-chuck set containing 2 twin-chucks and 1 central ring:
Note: These three parts are indissociable. In case of replacement, the three must be replaced simultaneously.

- Slide the chucks and the central ring on the pin and make sure that the flat ends of the chucks face their central ring.
- Slide the assembly into its lodging (axle).
- Screw the spindle onto the axle and tighten it using the wrench and the flat spanner.



2.4 Dismantling of the handpiece 49 E shell

- Unscrew and remove the nozzle. Pull the milled locking ring from the handpiece shell.

2.5 Dismantling of the handpiece 49 E shell

- Unscrew and remove the nozzle using the edge of the spanner, placed across the grooves in front of the axle.

FR

EN


DE

IT

ES

2.6 Reassembling the handpieces

- Slide the front of the wrench onto the rear of the axle by making sure that its two ins fit into the grooves of the axle's bearing.
- Whilst gently moving the milled locking ring on the handpiece shell, slide assembly into this latter.
- Push until the wrench touches the rear of the handpiece.
- Turn the wrench to its stop, then remove it.



III. ADVICE FOR USE

Always make sure that an instrument is in place before running the handpiece.


3.1 Insertion of the instrument

- Insert the instrument into the handpiece.
- Turn the milled locking ring counterclockwise.

3.2 Removal of the instrument

Open the chuck by fully turning the milled locking ring clockwise.

FR

EN

DE

IT

ES

3.3 Use of latch type drills

- Insert the depth pin assembly with its slit side forward into the handpiece.
- Remove the rear half of the depth pin assembly by unscrewing it.
- Your handpiece is ready for use with latch type drills.

IV. REPROCESSING OF THE HANDPIECE

(Follow instructions after every use.)

4.1 General information

(Follow instructions after every use.)

- The sterilization process does not affect the performance of dental handpieces and contra-angles. The end of the devices' service life is usually determined by damages due to use.
- Immediately after every use, withdraw excessive residues with a disinfecting wipe.
- Do not use products containing acetone, chlorine or bleach as disinfectants! If using disinfecting wipes or jets, follow the instructions of use for these products.
- It is recommended to reprocess the handpiece directly after every use. If this is not possible, immediately rinse the handpiece with water to prevent any impurities, blood or saliva from drying or coagulating inside the device.
- If blood has penetrated the handpiece, **immediately after use**, dismantle the handpiece and rinse it abundantly with cold water in order to eliminate the blood before it will have coagulated. Dry the handpiece, reassemble it and lubricate it with the "500 I" spray. **This process must be carried out before liquids like blood, physiological salt solution, saliva, etc. can dry or coagulate inside the handpiece.**
- Once a week, clean the dismantled parts using a degreasing liquid, dry them thoroughly before reassembling them and lubricate the assembled handpiece using a spray can "500 I".

4.2 Manual cleaning prior to sterilization

(Follow instructions after every use.)

The immersion of the handpiece during pre-disinfecting stage (immersion time must not exceed 15 minutes) is authorized in accordance with the following instructions:

- Use a detergent-disinfectant product specifically intended to the pre-disinfection of dental handpieces and contra-angles which does not contain soda based products or aldehyde.
- Comply with the concentration recommended by the pre-disinfectant manufacturer.
- Respect the immersion time recommended by the pre-disinfectant manufacturer.
- After immersion, rinse abundantly to eliminate the liquid / residues mixture contained in the handpiece.
- Immediately after rinsing, thoroughly dry the inside of the handpiece using compressed air.
- Before sterilization, lubricate the handpiece using an oil spray ("500 I" by MICRO-MEGA®).

The last three stages of this cleaning process (rinsing, drying and lubrication) may be performed by automatic devices available on the dental market, following the instructions of use of these devices.

MICRO-MEGA® recommends using X-Cidol® 2 as pre-disinfectant product.

FR

EN

DE

IT

ES

4.3 Automatic cleaning prior to sterilization


(Follow instructions after every use.)

In order to eliminate all handling, after each use on a patient place the handpiece in a device of the product line X-Cid®. These devices clean, disinfect and lubricate your handpieces and contra-angles without previous dismantling and guarantee a rapid bactericide, fungicide and virucide effect, ensuring the destruction of microbial agents such as HIV, HBV and Tuberculosis. The regular use of a device of the product line X-Cid® significantly extends the life span of your handpieces and contra-angles.

4.4 Maintenance (Lubrication)

(Follow instructions after every use.)

- Remove the instrument.
- Introduce the "500 I" spray into the rear of the handpiece.
- Press the spray a number of times lasting several seconds until the liquid coming out of the handpiece is clean.
- Make the handpiece turn in order to evacuate excess oil before use.
- During long treatments, repeat the lubrication process every ten minutes.

4.5 Sterilization

(Follow instructions after every use.)

- MICRO-MEGA® handpieces and contra-angles can be sterilized (wrapped in a sterilization bag) in a steam autoclave at 135°C (275°F) to 138°C (280.4°F) and 2.2 bar during 18 minutes.
- Perform sterilization in a sterilizer compliant with the standard EN 13060. Contact corrosion due to humidity may be possible.
- In case you want to sterilize several handpieces or contra-angles in one autoclave cycle, make sure that the autoclave's maximum load is not exceeded (cf. the autoclave's instructions of use).
- Immediately remove the device from the steam autoclave once the sterilization cycle is accomplished.
- Following manual or automatic preparation for sterilization and lubrication, place the handpiece in a hermetic sterilization bag (wrap the handpiece separately in a retractable sterilization bag):
 - Individual packaging: The sterilization bag must be of appropriate size for the handpiece and not be distended.
 - Batch processing: The handpiece may be placed in an instrument box provided to this effect or in a tray generally used for sterilization. Wrap up the trays using appropriate methods.
- Respect the sterilization bags' duration of storage indicated by the manufacturer.
- Quality and use of the wrapping material must be appropriate for sterilization and comply with the applicable standards.
- Open the sterilization bag just before use.

FR

EN

DE

IT

ES

V. REPAIRS

For service or repairs, send the complete handpiece.

Inspections and repairs may only be carried out by an approved dealer or by the manufacturer.

MICRO-MEGA® recommends annual service for all instruments.

VI. GUARANTEE AND AFTER SALES SERVICE

Your MICRO-MEGA® handpiece is an accurate medical device, produced under maximum quality and test conditions. To guarantee personal follow-up and traceability of its medical devices, MICRO-MEGA® would like you to return the guarantee card, duly filled in, within 10 days of your purchase. This return ensures that the procedure for after sales care of your handpiece will be instigated so that MICRO-MEGA® can give you full satisfaction with the service you expect and manage your guarantee correctly.

You can also register the guarantee for your handpiece on our website:

www.micro-mega.com

MICRO-MEGA® guarantees your handpiece for 12 months against any hidden manufacturing defect, excluding faults due to:

- Improper use.
- Inappropriate maintenance or application, not in accordance with our recommendations and instructions for use.
- Accidental damage (fall, impact, etc.).
- Attempted intervention or repair by personnel not approved by MICRO-MEGA®.
- Attempted product modification.

VII. SPECIFICATIONS

	49 E	49 EP	49 IS
Type	Assembly 1:1	Assembly 1:1	Assembly 1:1
Class (Directive CE 93/42)	IIA	IIA	IIA
Colour code	Blue	Blue	Blue
Maxi motor speed	40 000 rpm	40 000 rpm	40 000 rpm
Maxi instrument speed	40 000 rpm	40 000 rpm	40 000 rpm
Motor connection standard (coupling)	ISO Standard 3964	ISO Standard 3964	ISO Standard 3964
Instrument fixation	Mechanical locking	Mechanical locking	Mechanical locking
Instrument type in accordance with NF EN ISO 1797-1 (Instrument type)	Type 2	Type 2	Type 2
Instrument diameter (mm) in accordance with NF EN ISO 1797-1	2.35	2.35	2.35
Spray	Without spray	With external spray	With internal spray

FR

EN


DE

IT


ES

VIII. ACCESSORIES

8.1 References of accessories for straight handpiece PM 57 E


Depth pin assembly
Ref 59800197

Spray clip for 49 E
Ref 59806061

Spanner
Ref 90600013

Wrench
Ref 90600041

8.2 References of accessories recommended by MICRO-MEGA®

Ref 50700008

500 l Spray

supplied without nozzle

Ref 50700017

500 l Spray

supplied with 3 nozzles

Ref *Please contact us for further information.*
X-Cid® 2

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